
Sudan III, Practical grade

GRM991

Product Identifier

CAS No.	:	85-86-9
EC No.	:	201-638-4
C.I. No.	:	26100
Molecular Formula	:	C ₂₂ H ₁₆ N ₄ O
Molecular Weight	:	352.39
Synonym	:	Cerasin red
HS Code	:	3212 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Red to dark brown powder or solid
Solubility	:	20 ppm in toluene yields clear solution
Absorption maxima	:	506 - 512 nm (20 ppm in toluene)
Loss on drying (at 110°C, 2 hr)	:	<= 2.0%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	QK4250000
WGK	:	3