
Paraffin wax with ceresine**GRM9945B****Block packed in plastic container****Product Identifier**

CAS No.	:	8002-74-2
EC No.	:	232-315-6
HS Code	:	2712 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or white translucent blocks
Solubility	:	33.3 mg soluble in 1 mL of chloroform
Congeaing Point	:	60°C
Sulphated ash	:	<= 0.1%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	RV0350000
WGK	:	nwg