
Cedarwood oil, For Microscopy

GRM9971

Product Identifier

HS Code : 3301 29 16
Storage : Below 30°C
Shelf life : 2 years

Technical Specification

Appearance : Colourless to pale greenish yellow viscous liquid.
Solubility : 33.3 mg soluble in 1 mL of ethanol
Refractive index (n 20/D) : 1.513 - 1.524
Density (at 25°C) : 0.940 - 0.960 g/mL
Suitability for microscopy : Passes test

Safety Information

Hazard Pictogram(s) :



Signal Word : Warning
Hazard Statement(s) : H315
Precautionary Statement(s) : P280- P302 + P352- P321- P332 + P313- P362 + P364
UN No. : Not dangerous goods
Class : -
Packing Group : -
RTECS : FJ1520050
WGK : 1