



Technical Data

GW090

Concavity Slides Two Cavity

Product Number

Packing

GW090 : 1 PK

Product Information

Product Code : GW090
Product Name : Concavity Slides Two Cavity

Technical Specification

Description : Concavity slides having two polished spherical concavities are 16 mm in diameter x 0.5 mm deep, made of non corrosive glass with bevelled and polished edges and corners
Dimensions : 75 x 25 mm
Thickness : 1.4 mm

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516, Swastik Disha Business Park, Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India.
Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com