

## Blood Grouping Slide

**GW099**

Blood grouping slide 100 x 100 and 13 mm thick having 12 concavities each of 21 mm diameter. and 2.00 mm deep, surface of the slide frosted to provide a smooth marking, corners are bevelled and polished.

**Application :** Blood banks ,Medical and clinical laboratories etc.

Product Name	Product Code	Description	Size ( mm)	Thickness ( mm)	No. of Cavities
<b>Blood Grouping Slide</b>	<b>GW099</b>	Blood Grouping Slides are polished and frosted. Corners are bevelled.	100 X 100	13.00	12.00

### Product Features:

- Slides with frosted marking area.
- With polished and frosted
- Beveled edges.
- Packed in velvet lined box.
- Pack Size : 1 No.

### Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.