

---

## GW185 Homogeniser with plain pestle (P.P.) 20ml

---

Product Number	Packing
----------------	---------

GW185	: 1 No.
-------	---------

---

### Product Information

Product Code	: GW185
Product Name	: Homogeniser with plain pestle (P.P.) 20ml

---

### Technical Specification

Description	: Tissue homogenizer made up of borosilicate glass (grinded) vessel and plain pestle having PTFE head.
Capacity	: 20 mL
Use	: Used for homogenizing tissue.

---

### Transport Information

ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods