

HS001 HiSpark Alkaline Cleaning Solution, Biodegradable

Product Number Packing

HS001 : 500ML
HS001 : 5L

Product Information

Product Code : HS001
Product Name : HiSpark Alkaline Cleaning Solution, Biodegradable
HS Code : 3402 20 90
Other Information : Foam volume (5% in water) is not less than 70 ml
Shelf Life : 4 years

Technical Specification

Appearance : Colourless to pale yellow, clear to slightly turbid liquid.
Solubility : 33.3 mg soluble in 1 mL of water
pH (5% in water at 25°C) : 11.00 - 12.00
Density (at 20°C) : 1.065 - 1.075 g/mL

Risk and Safety Information

Storage Temperature(°C) : Store below 30°C

Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods