

## HS002 HiSpark Cleaning Solution, Neutral, Biodegradable

### Product Number

### Packing

HS002	: 500ML
HS002	: 5L

### Product Information

Product Code	: HS002
Product Name	: HiSpark Cleaning Solution, Neutral, Biodegradable
HS Code	: 3402 20 90
Other Information	: Foam volume (5% in water ) is not less than 80 ml
Shelf Life	: 4 years

### Technical Specification

Appearance	: Colourless to pale yellow, clear to slightly turbid liquid.
Solubility	: 33.3 mg soluble in 1 mL of water
pH (5% in water at 25°C)	: 6.00 - 8.00
Density (at 20°C)	: 1.057 - 1.067 g/mL

### Risk and Safety Information

Storage Temperature(°C)	: Store below 30°C
-------------------------	--------------------

### Transport Information

Marine Pollutant	: No
ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods