

HS003 HiSpark Cleaning Solution, Phosphate free, Biodegradable

Product Number

Packing

HS003	: 500ML
HS003	: 5L


Product Information

Product Code	: HS003
Product Name	: HiSpark Cleaning Solution, Phosphate free, Biodegradable
HS Code	: 3402 20 90
Shelf Life	: 4 years

Technical Specification

Appearance	: Colourless to pale yellow, clear to slightly turbid liquid
Solubility	: 33.3 mg soluble in 1 mL of water
Density (at 20°C)	: 1.058 - 1.066 g/mL
Free alkali (as NaOH)	: 1.5 - 2.5%
Phosphate (PO ₄)	: <= 0.002%

Risk and Safety Information

R-Phrase(s)	: 36/38
S-Phrase(s)	: 26
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	: 



Transport Information

Marine Pollutant	: No
ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods