

---

## Salt Testing Kit ( For Iodate )

**K052**

---

### Product Identifier

HS Code	:	3822 00 90
Storage	:	Below 30° C
Shelf life	:	1 Year
Packing	:	1PK (Sufficient for 1000 Salt samples)

---

### Technical Specification

Description	:	Kit box contents: Reagent Bottle (10 ml x 10 no.) Reagent is colourless non-viscous, translucent fluid.
Results	:	After addition of Reagent on salt sample colour changes are observed as:  Colour : Approximate Iodine content a)Light bluish violet colour : Less than 15 ppm b)Dark bluish violet colour : More than 15 ppm c)No change in colour : Absent

---

### Safety Information

UN No.	:	Not dangerous good
Class	:	-
Packing Group	:	-