
Salt Testing Kit (For Iodide)

K052A

Product Identifier

HS Code	:	3822 00 90
Storage	:	Below 30° C
Shelf life	:	1 Year
Packing	:	1PK (Sufficient for 500 Salt samples)

Technical Specification

Description	:	Kit box contents: Reagent K052A-1 (10 ml x 5 no.) and Reagent K052A-2 (10 ml x 5 no.) Reagents are colourless, non-viscous translucent fluid.
Results	:	After addition of Reagent on salt sample colour changes are observed as : Colour : Approximate Iodine content a) Light bluish violet colour : Less than 15 ppm b) Dark bluish violet colour : More than 15 ppm c) No change in colour : Absent

Safety Information

UN No.	:	Not dangerous good
Class	:	-
Packing Group	:	-