

---

## Soil Testing kit for Micronutrients

**K095S**

(Sufficient for 10 soil samples)

### Product Identifier

---

HS Code : 3822 00 90  
Storage : Below 30° C  
Shelf life : 1 Year

### Technical Specification

---

Description : The box contains the following

- 1) A booklet that describes the methods for determination of Copper (Cu), Zinc (Zn), Boron (B), Manganese (Mn), Iron (Fe) and Molybdenum (Mo) as well as interpretation of results.
- 2) Measuring cylinder (10 ml), Plastic beaker (50 ml), Glass rod, Funnel, Ordinary Filter Paper, Plastic soil measuring scoop (5g & 2g), Clean glass test tube, Small plastic spoon, stopper/rubber cork.
- 3) Reagent bottles - 17 Nos.
- 4) Colour chart

Results : The soil testing kit for Micronutrients is a box by which any one including the farmers can test the sample soil and obtain level of Copper (Cu), Zinc (Zn), Boron (B), Manganese (Mn), Iron (Fe) and Molybdenum (Mo) content in the soil.

### Safety Information

---

UN No. : Not dangerous good  
Class : -  
Packing Group : -