
Plant Tissue Test Kit (Sap Spot Test for B,Cu,Fe,Mn & Zn)

K107

(Sufficient for 50 soil samples)

Product Identifier

HS Code : 3822 00 90
Storage : Below 30° C
Shelf life : 1 Year

Technical Specification

Description : The box contains Compact flexible plastic carry case with PP buffer to hold all the reagent bottles and other items properly, Product insert that describes the methods for Plant Tissue Testing as well as interpretation of results, Reagent bottles(5 Nos),Treated Test Paper(3 Nos),Plastic cylinder (10 ml), Plastic dropper,Glass Test tube,5 ml(1 Nos),Pricking pin, plastic bag plain, & perforated etc

Results : The Plant Tissue test kit is a box by which any one including the farmers can test the Plant Tissue (leaf) sample for B,Cu, Fe,Mn & Zn values.

Safety Information

UN No. : Not dangerous good
Class : -
Packing Group : -