

## Stainless Steel caps

LA003, LA004,  
LA005, LA735.

This cap is made of stainless steel and having firm grip for rimless glass test tube.

**Application :** In Microbiology laboratory, Culture preparation and Various laboratories.

Product Name	Product Code	Size Approx. OD mm
Firm Grip S.S. Cap for rimless glass test tube.	LA735	15.00
Firm Grip S.S. Cap for rimless glass test tube.	LA003	18.00
Firm Grip S.S. Cap for rimless glass test tube.	LA004	25.00
Firm Grip S.S. Cap for rimless glass test tube.	LA005	38.00

### Product Features:

- Firm Grip.
- Pack Size : 1X 50 No / pack

### Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.