

Product Info

1. TabSpin®006 is compact and cost effective, allowing any workstation to be equipped with a “personal” mini centrifuge.
2. It is supplied with one rotor suitable for 6-1.5/2.0ml tubes, 6-0.5ml tubes and 2-8x0.2ml PCR tubes, and two speed options (4,000 & 6,000rpm).
3. The electronic braking can reach quick deceleration to reduce handling time.



TabSpin®006 LA1089

Features

- 3-in-1 Rotor
- Compact and cost effective
- Two speed options: 4,000 & 6,000rpm
- TabSpin®006 is supplied with 1 rotor: 6 x 1.5/2.0ml, 6 x 0.5ml tubes and 2-8 x 0.2ml PCR Strip rotor

Technical Data

| Model | TabSpin®006 |
|--------------------|--|
| Capacity | 6 x 1.5/2.0ml, 6 x 0.5ml tubes and 2-8 x 0.2ml PCR Strip rotor |
| Two Speed options | 4,000rpm & 6,000rpm |
| Maximum RCF | 1,200 & 2,000 x g |
| Weight | 0.7kg |
| Electrical | AC100~240V, 50/60Hz |
| Dimensions (DxWxH) | 135 x 127 x 120mm |

