

LA1126 : Micropipette Tips- 200µl, Yellow, Nuclease free

Product Number : Packing

LA1126 : 1X1000NO

LA1126 : 5X1000NO

Product Information

Product Code : LA1126

Product Name : Micropipette Tips- 200µl, Yellow, Nuclease free

Technical Specification

Description : Pre-sterile universal tips. Tips are Nuclease free.

Clarity : Clear transparent without any particles and Fibers

Capacity : 200 µL

Material : Medical Grade Polypropylene

Color : Yellow

Graduations : No

Autoclavable : Yes

Height : Approx. 50.6 mm

Inner Diameter : Approx. 5.2 mm

Outer Diameter : Approx. 7.4 mm

Product alignment : Straight without any bending.

Use : Used to dispense liquid samples.

Compatible HiMedia's

Pipette codes : LA599, LA600, LA612, LA613, LA615, LA618, LA959,

LA965, LA966, LA967, LA958


Please refer disclaimer Overleaf.

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

Technical Assistance

At HiMedia we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at mb@himedialabs.com.



Storage temperature



Do not use if package is damaged



HiMedia Laboratories Private Limited,
Reg. Off: Plot No. C-40, Road No. 21Y,
MIDC, Wagle Industrial Area, Thane,
(West) 400604, Maharashtra, INDIA.
Web: www.himedialabs.com



06/2025

PILA1126_0/0622

LA1126-00

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: techhelp@himedialabs.com Website: www.himedialabs.com