

**LA1128 : Micropipette tips-200µl, clear, Nuclease free**

**Product Number : Packing**

LA1128 : 1X1000NO

**Product Information**

Product Code : LA1128

Product Name : Micropipette tips-200µl, clear, Nuclease free

**Technical Specification**

Description : Pre-sterile universal tips. Tips are Nuclease free.

Clarity : Clear transparent without any particles and Fibers

Capacity : 200 µL

Material : Medical Grade Polypropylene

Color : clear

Graduations : No

Autoclavable : Yes

Height : Approx. 50.2 mm

Inner Diameter : Approx. 5.3 mm

Outer Diameter : Approx. 7.2 mm

Autoclavable : Yes

Product alignment : Straight without any bending.

Use : Used to dispense liquid samples.

Compatible HiMedia's

Pipette codes : LA599, LA600, LA612, LA613, LA615, LA618, LA959,  
LA965, LA966, LA967, LA958

**Please refer disclaimer Overleaf.**

## Transport Information


---

ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods

## Technical Assistance

---

At HiMedia we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at [mb@himedialabs.com](mailto:mb@himedialabs.com).



Storage temperature



Do not use if package is damaged



HiMedia Laboratories Private Limited,  
Reg. Off: Plot No. C-40, Road No. 21Y,  
MIDC, Wagle Industrial Area, Thane,  
(West) 400604, Maharashtra, INDIA.  
Web: [www.himedialabs.com](http://www.himedialabs.com)



11/2025

PILA1128\_0/1122

LA1128-00

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.  
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: [techhelp@himedialabs.com](mailto:techhelp@himedialabs.com) Website: [www.himedialabs.com](http://www.himedialabs.com)