

LA1128A : Micropipette tips-200µl, clear, Nuclease free

Product Number : Packing

LA1128A : 1X1000NO

LA1128A : 5X1000NO

Product Information

Product Code : LA1128A

Product Name : Micropipette tips-200µl, clear, Nuclease free

Technical Specification

Description : Non sterile universal tips. Tips are Nuclease free.

Clarity : Clear transparent without any particles and fibers

Capacity : 200 µL

Low retention : No

Material : Medical Grade Polypropylene

Color : clear

Graduations : Yes

Autoclavable : Yes

Sterile : No

Height : Approx. 58 mm

Inner Diameter : Approx. 5.2 mm

Outer Diameter : Approx. 7.3 mm

Autoclavable : Yes

Product alignment : Straight without any bending.

Use : Used to dispense liquid samples.

Please refer disclaimer Overleaf.

Compatible HiMedia's

Pipette codes : LA599, LA600, LA612, LA613, LA615, LA618, LA959,
LA965, LA966, LA967, LA958

Transport Information

ADR/RID : Not Dangerous Goods









IMDG : Not Dangerous Goods

IATA : Not Dangerous Goods

Technical Assistance

At HiMedia we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at mb@himedialabs.com.

Symbols:

	Manufacturer		Do not use if package is damaged
	Batch code		Temperature limit
	Date of manufacture (YYYY-MM)		Consult instructions for use
	Use-by date (YYYY-MM)		Catalogue number

Identification No.: PILA1128A

Rev No.: 00

Date of Issue.: 2025-05

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: techhelp@himedialabs.com Website: www.himedialabs.com