
LA213 S.S. Test Tube Racks 18 - 5 x 10 x90mm ht.

Product Number

Packing

LA213 : 1 NO

Product Information

Product Code : LA213
Product Name : S.S. Test Tube Racks 18 - 5 x 10 x90mm ht.

Technical Specification

Description : Stainless steel test tube racks made out of 3 mm diameter wire
Dimensions : a) Height : 90 mm b) Diameter : 18 mm
Capacity : 5 X 10
Use : To hold test tubes

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.