

LA222

S.S.Test Tube Racks 38-5x4x100mm ht. Sh



Product Number**Packing**

LA222 : 1NO

Product Information

Product Code : LA222
Product Name : S.S.Test Tube Racks 38-5x4x100mm ht. Sh

Technical Specification

Description : Stainless steel test tube racks made out of 20g sheet.
Dimensions : a) Height : 100 mm b) Diameter : 38 mm
Capacity : 5 X 4
Use : To hold test tubes

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods