



# Technical Data

**LA225**

**S.S. Test Tube Racks 38;5 x 4 x100mm ht.**

---

## Product Number

## Packing

LA225 : 1 NO

---

## Product Information

Product Code : LA225  
Product Name : S.S. Test Tube Racks 38;5 x 4 x100mm ht.

---

## Technical Specification

Description : Stainless steel test tube racks made out of 3 mm diameter wire.  
Diameter of holding test tubes : 38 mm  
Height : 100 mm  
Diameter of holding test tubes : 5x4  
Use : To hold test tubes

---

## Transport Information

ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

---

HiMedia Laboratories Pvt. Ltd. A-516,Swastik Disha Business Park,Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India.  
Customer care No.: 022-6147 1919 Email: [techhelp@himedialabs.com](mailto:techhelp@himedialabs.com) Website: [www.himedialabs.com](http://www.himedialabs.com)