

## S.S. Test Tube Basket

LA277

It can withstand alkali and acid solutions. Made out of dia 3.00 mm wire and available in square shape. This basket is used to washing test tubes and has high strength corrosion.

**Size :** 6" X 6" X 6" .(L X W X H)

**Application :** Hospitals , Microbiology laboratory, Molecular Biology laboratory and various other laboratories.

**Pack Size :** 1 No. / Pack.

### Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.