



## Auto HiPette™

**Product code: LA692, LA692A, LA692B**

### Product description:

Auto HiPette™ provides a convenient aid for dispensing liquid from glass and serological pipettes.

### Features:

- Has a safety valve and hydrophobic filters that provide protection against fluid penetration.
- Three different speed modes- High, Low and Gravity Drain Provide better control of speed (convenience of fine-tuning pump speed).
- Cadmium-free environment-friendly NiMH batteries (extended operational life).
- Ease of changing batteries.
- Has a low battery indicator

### Specifications:

Product	Code	Packing
Auto HiPette™	LA692-1No	1 No
Auto HiPette™ Battery	LA692A-2No	2 No
Accessory for Auto HiPette™	LA692B-1PK	1 PK

### Disclaimer :

Revision No.: 02/2022

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.