

Stainless Steel Scalpel Holder No.-3

LA723

New scalpel holder made with stainless steel grade SS 410. With excellent expertise, we supply a supreme quality range of Scalpel Blades that are made from stainless steel. They are used in medical sector and in different surgeries.

Applications : Hospitals, Health care centers, Nursing homes, Medical laboratory, Microbiology laboratory , Plant tissue culture work and various laboratories.

Product Features:

- Touch sharpness
- Light weight with long neck
- Free from burrs
- Fine finishes
- Pack Size : 2 Nos

Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.