

LA749XL : Barrier Tips XL, Max Capacity 10µl

Product Number : Packing

LA749XL : 1X96NO

LA749XL : 10X96NO

Product Information

Product Code : LA749XL

Product Name : Barrier Tips XL, Max Capacity 10µl

Technical Specification

Description : Nuclease free/ Pre-sterile universal low retention barrier tips.

Clarity : Clear transparent without any particles and Fibers

Capacity : 10 µL

Material : Tips – Medical Grade Polypropylene
Barrier – Medical Grade HDPE

Color : Clear, natural

Graduations : Yes

Autoclavable : NO

Height : Approx. 44.7 mm

Inner Diameter : Approx. 3.4 mm

Outer Diameter : Approx. 5.8 mm

Product alignment : Straight without any bending.

Use : Used to dispense liquid samples.

Compatible HiMedia's

Pipette codes : LA611, LA617, LA955, LA964


Please refer disclaimer Overleaf.

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

Technical Assistance

At HiMedia we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at mb@himedialabs.com.



Storage temperature



Do not use if package is damaged



HiMedia Laboratories Private Limited,
Reg. Off: Plot No. C-40, Road No. 21Y,
MIDC, Wagle Industrial Area, Thane,
(West) 400604, Maharashtra, INDIA.
Web: www.himedialabs.com



06/2025

PILA749XL_0/0622

LA749XL-00

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: techhelp@himedialabs.com Website: www.himedialabs.com