

LA792 Assay Cylinder

Product Number	Packing
LA792	: 1 x 10 No.
LA792	: 1 x 20 No.

Product Information

Product Code	: LA792
Product Name	: Assay Cylinder

Technical Specification

Appearance	: Assay cylinders are stainless steel cylinders.
Dimension	: Outer diameter : 8mm ± 0.1mm Inner diameter : 6mm ± 0.1mm Length: 10mm ± 0.1mm

Transport Information

ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods