

LA859C : Barrier Tips, 1000µl, Max capacity 1250µl

Product Number : Packing

LA859C : 1X96NO

LA859C : 10X96NO

Product Information

Product Code : LA859C

Product Name : Barrier Tips, 1000µl, Max capacity 1250µl

Technical Specification

Description : Pre-sterile universal low retention barrier tips.

Tips are Nuclease free.

Clarity : Clear transparent without any particles and Fibers

Capacity : 1250 µL

Material : Tips – Medical Grade Polypropylene

Barrier – Medical Grade HDPE

Color : Clear, natural

Graduations : Yes

Height : Approx. 101.4 mm

Inner Diameter : Approx. 7.3 mm

Outer Diameter : Approx. 8.8 mm

Product alignment : Straight without any bending.

Use : Used to dispense liquid samples.

Compatible HiMedia's

Pipette codes : LA614, LA616, LA960, LA606, LA607, LA608

Please refer disclaimer Overleaf.









Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

Technical Assistance

At HiMedia we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at mb@himedialabs.com.

Symbols:

| | | | |
|---|----------------------------------|--|----------------------------------|
|  | Manufacturer |  | Do not use if package is damaged |
|  | Batch code |  | Temperature limit |
|  | Date of manufacture (YYYY-MM) |  | Consult instructions for use |
|  | Use-by date (YYYY-MM) |  | Catalogue number |

Identification No.: PILA859C

Rev No.: 01

Date of Issue.: 2025-07

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: techhelp@himedialabs.com Website: www.himedialabs.com