

Stainless Steel Petri Dish Carrier

LA276, LA774 & LA860

Stainless steel petri plate dish carriers to facilitate the handling and the transport of petri plates around the laboratory and on site.

Application : Useful for safe handling , transport and incubation of Petri dishes around the laboratory .

Product Name	Product Code	Description	Size (mm)
S.S Petri Dish Carrier	LA276	Stainless steel Petri dish carrier with handle, made out of S.S. 304 .	For dia of 100 mm dish.
S.S Petri Dish Carrier	LA774	Stainless Steel Petri dish carrier without handle, made out of S.S. 304.	For Dia 55.00 mm Petri plate
S.S Petri Dish Carrier	LA860	Stainless Steel Petri dish carrier with handle made out of S.S.304.	For Dia 55.00 mm Petri plates

Product Features

- The stainless steel construction, no corners to clean.
- Swinging handle for easy manipulation is what the operators are asking for.
- Used for sterilization and storage.
- Pack Size : 1 No./ Pack

Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.