

S.S. Spirit Lamp

LA275, LA873 & LA875.

Spirit Lamp with cap made from Stainless steel supplied with wick.

Application : Laboratory Use.

Product Name	Product Code	Description	Size (mm) (Ht X Dia)	Weight (g)	Capacity (ml)
S.S. Spirit Lamp	LA275	Stainless steel spirit lamp with cap and wick.	96.60 X 51.00	147.00	95.00 Spirit.
S.S. Spirit Lamp	LA873	Stainless steel spirit lamp with cap and wick is light in weight and made out of S.S 304.	63.00 X 63.00	55.50	125.00 Spirit
S.S. Spirit Lamp	LA875	Stainless steel spirit lamp with cap and wick is light in weight and made out of economy grade.	75.50 X 61.00	55.00	170.00 Spirit

Pack size : 1 No /Pack.

Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.