
LA876 Triangular S.S. Spreader with Bakelite Handle

Product Number	Packing
LA876	: 1 X 10 No.

Product Information

Product Code	: LA876
Product Name	: Triangular S.S. Spreader with Bakelite Handle

Technical Specification

Description	: Triangular stainless steel spreader has a triangular end with bakelite handle featuring smooth surfaces throughout, with a straight spreading bar.
Use	: Easily perform inoculation of petri dishes.

Transport Information

ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods