



# Product Information

**LA878**

**Varivol II Micropipette - 10000**

---

**Product Number**

**Packing**

LA878 : 1 No.

---

**Product Information**

Product Code : LA878  
Product Name : Varivol II Micropipette - 10000

---

**Technical Specification**

Description : Green coloured Variable volume microlitre pipettes wherein volumes can be continuously adjusted by rotating the plunger or thumb button. It has click-stop digital system, reverse pipetting technique and separate, streamlined tip ejector. Each pipette has been tested, calibrated and comes with its own calibration certificate.

Capacity : 1000 - 10,000  $\mu$ L

Use : Used in laboratory to dispense liquid samples.

---

**Transport Information**

ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919  
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

---

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein