

## Petri Dish Stand For Environmental Monitoring

LA889 & LA890.

Ideal for exposing plates in Clean Rooms for bacterial count.



Product	Product Code	Description	Size (cm)
<b>Petri Dish Stand For Environmental Monitoring- Parallel</b>	LA889	Ideal for exposing plates in clean rooms for bacterial counts.	20 x 11 x 30
<b>Petri Dish Stand For Environmental Monitoring- Inclined</b>	LA890	Ideal for exposing plates in clean rooms for bacterial counts.	20 x 11 x 30/35

### Product Features :

- Full SS (304) construction.
- Stackable space saving design.
- Ideal for storage when not in use.
- Ideal height of app 30 cm - to place below risers which are at 45 cm height.

### Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.