



Product Information

LA896 S.S. Petri Dish Carrier 55 mm, with handle

Product Number**Packing**

LA896 : 1 No.

Product Information

Product Code : LA896
Product Name : S.S. Petri Dish Carrier 55 mm, with handle

Technical Specification

Description : Stainless steel Petri dish carrier made out of SS 304.
Size of Petri Plate : For 55 mm diameter
Capacity to hold Petri Plate : For 11 plates
Temperature tolerance test : Withstands steam sterilization at 121°C for 30 minutes and not to be used in dry heat sterilization.
Use : For laboratory use as Petri dish carrier.

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods