

LA937B : Freeze Tag, Blue Dots

Product Number : Packing

LA937B : 1X1000NO

LA937B : 10X1000NO

Product Information

Product Code : LA937B

Product Name : Freeze Tag, Blue Dots

Technical Specification

Description : Freeze Tags are made from a glossy high quality polypropylene film with a highly aggressive permanent rubber-based adhesive used to label vials, Cryo Tubes, microtubes, microplates, Freezer boxes and Freezer racks

Material : Polypropylene

Sheet quality : 45 GSM

Color : Blue

Diameter : 9.5mm

Printer friendly : Yes

Temperature Tolerance

Range : -40°C to +70°C

Use : Used to label difficult substrates with rough and textured surfaces.


Please refer disclaimer Overleaf.

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

Technical Assistance

At HiMedia we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at mb@himedialabs.com.



Storage temperature



Do not use if package is damaged



HiMedia Laboratories Private Limited,
Reg. Off: Plot No. C-40, Road No. 21Y,
MIDC, Wagle Industrial Area, Thane,
(West) 400604, Maharashtra, INDIA.
Web: www.himedialabs.com



07/2025

PILA937B_1/0722

LA937B-01

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: techhelp@himedialabs.com Website: www.himedialabs.com