

**LA940W : Freeze Tag, White Dots**

**Product Number : Packing**

---

LA940W : 1X1000NO

LA940W : 10X1000NO

**Product Information**

---

Product Code : LA940W

Product Name : Freeze Tag, White Dots

**Technical Specification**

---

Description : Freeze Tags are made from a matt opaque synthetic paper material with a specially formulated clear permanent pure acrylic adhesive featuring standard initial tack and bond strength to a wide range of substrates even at extreme high temperatures. These tags are chemically inert and resist to most organic solvents and caustic agents. They adhere to most plastics, glass and metals without cracking, peeling or degrading.

Material : Polypropylene

Quantity : 1000 labels/ pack, Roll form

Sheet quality : 66 GSM

Color : White

Diameter : 9.5 mm

Printer friendly : Yes

Temperature Tolerance

Range : -196°C to +121°C

Use : Used to label difficult substrates with rough and textured surfaces.

---

**Please refer disclaimer Overleaf.**

## Transport Information

---

ADR/RID : Not Dangerous Goods

IMDG : Not Dangerous Goods









IATA : Not Dangerous Goods

## Technical Assistance

---

At HiMedia we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at [mb@himedialabs.com](mailto:mb@himedialabs.com).

### Symbols:

	Manufacturer		Do not use if package is damaged
	Batch code		Temperature limit
	Date of manufacture (YYYY-MM)		Consult instructions for use
	Use-by date (YYYY-MM)		Catalogue number

Identification No.: PILA940W

Rev No.: 01

Date of Issue.: 2025-03

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.  
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: [techhelp@himedialabs.com](mailto:techhelp@himedialabs.com) Website: [www.himedialabs.com](http://www.himedialabs.com)