



Soyabean Casein Digest Medium w/BCP

LQ148XX

Intended Use:

Recommended for cultivation of wide variety of microorganisms with the addition of carbohydrates it can be used for fermentation studies.

Composition**

Ingredients	g / L
Tryptone	17.000
Soya peptone	3.000
Sodium chloride	5.000
Dipotassium hydrogen phosphate	2.500
Dextrose (Glucose)	2.500
Bromocresol purple	0.010
Final pH (at 25°C)	7.3±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Label the ready to use LQ148XX bottle. Inoculate 50-100 cfu sample and Incubate at specified temperature and time.

Principle And Interpretation

Soyabean Casein Digest Medium is a nutritious medium that will support the growth of a wide variety of microorganisms, including common aerobic, facultative and anaerobic bacteria and fungi (1-4). It can also be used as a general, all purpose cultivation medium (3). When tested for the growth of organisms in presence of indicator like bromocresol purple, the colour of the medium changes from purple to yellow. With the addition of carbohydrates it can be also used for the fermentation studies of fastidious and non-fastidious organisms.

Tryptone and soya peptone provides necessary amino acids and other complex nitrogenous substances. Dextrose serves as an energy source. Sodium chloride maintains the osmotic equilibrium. Dipotassium phosphate acts as a buffer to control pH.

Type of specimen

Environmental monitoring- Sterility testing in manufacturing process.

Specimen Collection and Handling

For environmental monitoring, follow appropriate techniques for sample processing in case of viscous materials as mentioned under sterility (1,4).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

1. Further biochemical testing is required on colonies of pure culture for complete identification.
2. Individual organisms differ in their growth requirement and may show variable growth patterns on the medium.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Sterile clear Soyabean Casein Digest Medium w/BCP in glass bottle.

Please refer disclaimer Overleaf.

Colour

Purple coloured clear solution

Quantity of Medium

20 ml of medium in glass bottle.

pH

7.10-7.50

Sterility Check

Passes release criteria

Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 18-48 hours (*Candida albicans* incubated at 20-25°C for 2-7 days).

Organism	Inoculum (CFU)	Growth	Colour of medium
<i>Candida albicans</i> ATCC 10231 (00054*)	50-100	luxuriant	yellow
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> ATCC 25923 (00034*)	50-100	good-luxuriant	yellow
<i>Escherichia coli</i> ATCC 25922 (00013*)	50-100	good-luxuriant	yellow
** <i>Bacillus spizizenii</i> ATCC 6633 (00003*)	50-100	good-luxuriant	yellow

Key : (*) Corresponding WDCM numbers.

**Formerly known as *Bacillus subtilis* subsp. *spizizenii*

Storage and Shelf Life

On receipt store between 15-30°C. Use before expiry date on the label. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (5,6).

Reference

- Forbes, Sahn and Weissfeld. 1998. Bailey & Scotts diagnostic microbiology, 10th ed. Mosby, Inc. St. Louis, Mo.
- Fredette and Forget. 1961. The sensitivity of several media to small inocula. Extract from a paper presented at the Canadian Society of Microbiology Annual Meeting, June 12-15. Kingston, Ontario, Canada.
- MacFaddin. 1985. Media for isolation-cultivation-identification-maintenance of medical bacteria, vol. 1. Williams & Wilkins, Baltimore, Md.
- Marshall (ed.). 1993. Standard methods for the examination of dairy products, 16th ed. American Public Health Association, Washington, D.C.
- Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
- Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W.(2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Revision : 00/2025

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.