

High Salt Peptone Yeast Extract Agar

Intended Use:

Recommended for isolation, cultivation and confirmation of salt tolerant *Vibrio* species in accordance with ISO Committee under specifications ISO/DIS 8914 : 1990.

Composition**

Ingredients	Gms / Litre
Peptone	10.000
HM extract#	2.000
Yeast extract	6.000
Sodium chloride	30.000
L-Cysteine hydrochloride	0.300
Glucose (Dextrose)	2.000
Agar	15.000
Final pH (at 25°C)	7.5±0.2
**Formula adjusted, standardized to suit performance parameters	

Equivalent to Meat extract

Directions

Suspend 65.3 grams in 1000 ml purified/distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri Plates.

Principle And Interpretation

Vibrios are fairly easy to isolate from both clinical and environmental materials, though some species may require growth factors and /or vitamins. *Vibrio parahaemolyticus* is the leading cause of bacterial diarrhea associated with the consumption of contaminated food products. Media can be made selective for *Vibrios* by adding appropriate selective agents (1). High concentrations of NaCl and alkaline pH have also been used to select certain *Vibrio* species, based on their ability to grow at pH values above 8.0 and at 3% or higher concentrations of NaCl.

Vibrio cholerae is a non-halophilic *Vibrio*, which cannot grow in media with a concentration of sodium chloride greater than 5-6% and is able to grow in media lacking NaCl (2). High Salt Nutrient Agar/High Salt Peptone Yeast Extract Agar is recommended for the isolation, cultivation and confirmation of salt-tolerant *Vibrio* species in products intended for human consumption or animal feeding stuffs in accordance with ISO Committee under specification ISO/DIS 8914:1990 (3).

HM extract, yeast extract, L-cysteine hydrochloride and peptone are sources of carbon, nitrogen, vitamins and minerals. Glucose is the fermentable carbohydrate. Sodium chloride maintains the osmotic balance of the medium and provides the essential ions.

Type of specimen

Clinical samples -Faeces

Specimen Collection and Handling

Inoculate 25 grams of test portion to 225 ml Salt Polymyxin Broth Base (M821I) and 225 ml Alkaline Peptone Water (M618I). Incubate the two broths at 35-37°C for 7-8 hours. After incubation, inoculate a loopful from M821I onto TCBS Agar (M189), Tryptone Sucrose Tetrazolium Agar Base (M1217) and High Salt Nutrient Agar (M1218). Repeat the plating procedure for M618I. Incubate the plates at 35-37°C for 20-24 hours. Confirm presumptive *Vibrio* colonies by performing the biochemical tests. This can be performed by inoculation into High Salt Peptone Yeast Extract Agar (M1219). This medium can be used to differentiate between aerobic and anaerobic growth. After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions :

In Vitro diagnostic Use. For professional use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

M1219

Limitations :

1. Individual organisms differ in their growth requirement and may show variable growth patterns on the medium.

2. Each lot of the medium has been tested for the organisms specified on the COA. It is recommended to users to validate the medium for any specific microorganism other than mentioned in the COA based on the user's unique requirement.

3. Further biochemical and serological tests must be carried out for further identification.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Light yellow coloured, clear to slightly opalescent gel forms in tubes

Reaction

Reaction of 6.53% w/v aqueous solution at 25°C. pH : 7.5±0.2

pН

7.30-7.70

Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organism	Inoculum (CFU)	Growth
<i>Vibrio cholerae</i> ATCC 15748	50-100	good-luxuriant
<i>Vibrio parahaemolyticus</i> ATCC 17802 (00037*)	50-100	good-luxuriant

Key : *Corresponding WDCM numbers.

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 20-30°C.Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition.Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (4,5).

Reference

1. Collee J. G., Fraser A. G., Marmion B. P., Simmons A., (Eds.), Mackie and McCartney Practical Medical

Microbiology, 1996, 14th Edition, Churchill Livingstone.

2.Bruno Gomez-Gil and Roque A., Isolation, Enumeration and Preservation of the Vibrionaceae, Thompson F. L., Austin B. and Swings J., The Biology of Vibrios, ASM press.

3. International Organization for Standardization (ISO), 1990, Draft ISO/DIS 8914:1990.

4. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.

5.Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Revision :03/2022

IVD	In vitro diagnostic medical device
(€	CE Marking
-30°C	Storage temperature
	Do not use if package is damaged
	HiMedia Laboratories Pvt. Limited, Plot No.C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) -400604, MS, India
EC REP	CE Partner 4U ,Esdoornlaan 13, 3951 DB Maarn The Netherlands, <u>www.cepartner</u> 4u.eu

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Corporate Office : Plot No.C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) - 400604, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com