



Technical Data

Yersinia Identification Broth Base

M1221

Intended use

Yersinia Identification Broth Base with addition of Urea is recommended for identification of *Yersinia*.

Composition**

Ingredients	Gms / Litre
L-Tryptophan	3.000
Potassium dihydrogen phosphate	1.000
Dipotassium hydrogen phosphate	1.000
Sodium chloride	5.000
Phenol red	0.025
Final pH (at 25°C)	6.9±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 10.02 grams in 1000 ml purified / distilled water. Heat if necessary to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C and aseptically add Urea solution (FD048). Mix well and dispense in sterile tubes or flasks as desired.

Principle And Interpretation

There are three species of *Yersinia* with unquestionable pathogenicity for humans *Yersinia pestis*, *Yersinia pseudotuberculosis* and *Yersinia enterocolitica*. Among these, *Y. enterocolitica* is usually associated with foodborne gastroenteritis. It is the most common species of *Yersinia* recovered from clinical specimens. The portal of entry in humans is the oral digestive route, with infection occurring in the terminal ileum (6). Yersinia Identification Broth Base is recommended for the identification of *Yersinia* from food and animal feeds by the ISO Committee (3).

L-Tryptophan serves as a base to test indole reaction. Phosphates buffer the medium while sodium chloride maintains the osmotic equilibrium of the medium. Phenol red is the pH indicator dye. Urea (FD048) is broken down by enzyme urease to yield ammonia. Ammonia increases the pH of the medium towards alkalinity, consequently making the phenol red indicator dye to change from an orange-red to a pink-violet colour.

Inoculate the test sample in PSB Broth (M941) and ITC Broth (M1220) for enrichment. After incubation at 25°C for 2-3 days, inoculate onto Yersinia Selective Agar Base (M843). Presumptive Yersinia colonies are confirmed biochemically by inoculating into Yersinia Identification Broth Base (M1221) (3).

Type of specimen

Animal feed and feeding stuff

Specimen Collection and Handling:

For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (2,5,6). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions :

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations :

1. Further isolation and biochemical tests must be carried out for confirmation

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Please refer disclaimer Overleaf.

Quality Control

Appearance

Light yellow to light pink homogeneous free flowing powder

Colour and Clarity of prepared medium

Orange-red coloured clear solution without any precipitate.

Reaction

Reaction of 1.0% w/v aqueous solution at 25°C. pH : 6.9±0.2

pH

6.70-7.10

Cultural Response

Cultural characteristics observed, after an incubation at 30-32°C for 18-24 hours with added Urea solution (FD048).

Organism	Inoculum (CFU)	Growth	Urease production	Colour change of medium
<i>Yersinia enterocolitica</i> ATCC 27729	50-100	luxuriant	positive reaction	Orange-red to cerise

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (4, 5).

Reference

1. American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington D.C.
2. Downes F. P. and Ito K., (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 4th Ed., American Public Health Association, Washington, D.C.
3. International organization for standardization, (ISO), 1994, Draft ISO 10273. Microbiology of food and animal feeding stuffs--Horizontal method for the detection of presumptive pathogenic *Yersinia enterocolitica*.
4. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
5. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W.(2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
6. Koneman E. W., Allen S. D., Janda W. M., Schreckenberger P. C., Winn W. C. Jr., 1992, Colour Atlas and Textbook of Diagnostic Microbiology, 4th Ed., J. B. Lippincott Company.

Revision : 03/2018

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.