



Technical Data

Bryant and Burkey Medium

M1385

Intended Use:

Recommended for detecting and enumerating spores of lactate fermenting *Clostridium* in milk and dairy products.

Composition**

Ingredients	Gms / Litre
Tryptone	15.000
Yeast extract	5.000
HM peptone B #	7.500
Sodium acetate	5.000
L-Cystine dihydrochloride	0.500
Resazurin	0.0025
Final pH (at 25°C)	5.9±0.2

**Formula adjusted, standardized to suit performance parameters

* Equivalent to Beef extract

Directions

Suspend 33.0 grams in 1000 ml purified/distilled water containing 5 grams of Sodium lactate. Heat if necessary to dissolve the medium completely. Mix well and dispense into tubes or flasks as desired. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

Principle And Interpretation

Bryant and Burkey Medium is used to enumerate the spores of gasogenic clostridia that are responsible for swelling of cheese in dairy industry (1, 2). In normal conditions of use, the medium allows the growth of other microorganisms also which are not directly related to cheese alteration, e.g. *Clostridium sporogenes* or *Clostridium butyricum*. Germination and growth of butyric acid bacteria (BAB) causes severe defects in cheese with silage being the main source of BAB spores in cheese milk (6). Clostridia spores are heat resistant and therefore can contaminate cheese brines. The gas produced by the growth of clostridia swells the cheese and is responsible for defect known as butyric swelling, resulting in bad taste. The main species causing this butyric swelling defect is *Clostridium tyrobutyrium*.

Recommended technique for estimation of *Clostridium* is to enumerate the spores by the MPN method. Resazurin is a redox indicator and monitors the oxygen level turning from pink under aerobic conditions to colourless under anaerobic conditions. The nutrient composition of the basal medium, particularly of tryptone, yeast extract, HM peptone B and L-cysteine HCl help in rapid growth of *Clostridium* species. Sodium acetate promotes spore germination (5) and improves the selectivity of the medium. Sodium lactate is fermented under anaerobic conditions by *Clostridium tyrobutyricum*. A pink colour indicates the presence of oxygen in the media which turns colourless on boiling.

Type of specimen

Dairy samples

Specimen Collection and Handling

Test sample must be previously decontaminated by heating up for 10 minutes at 75°C in order to destroy all the vegetative forms while only leaving the spores alive. The tubes of the medium are then boiled to regenerate anaerobiosis and cooled down to 25-30°C. These tubes are inoculated with the sample and overlaid with 2 cm of sterile paraffin and incubated for upto 7 days at 37°C. Examine the tubes after every 48 hours. A raised paraffin plug indicates gas formation. Tubes will be declared positive if they show clear gas production and MPN index is used to calculate the number of Clostridia. After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions :

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations :

1. Proper anaerobic conditions to be provided for the growth of Clostridia.
2. Some strains may show poor growth due to nutritional variations.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Colour and Clarity of prepared medium

Light amber coloured, clear solution when hot, becomes pink-red upon cooling

Reaction

Reaction of 3.3% w/v aqueous solution at 25°C. pH : 5.9±0.2

pH

5.70-6.10

Cultural Response

Cultural characteristics observed under anaerobic conditions after an incubation at 35-37°C for 6 days.

Organism	Growth	Gas production
----------	--------	----------------

<i>Clostridium tyrobutyricum</i> ATCC 25755	good	positive
--	------	----------

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 15-25°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (3,4).

Reference

1. Bergeres J. L. and Sivela S., 1989, Detection and enumeration of Clostridial spores related to cheese quality-Classical and new methods, FIL-IDF Bull. 51:18-23.
2. Bryant M. P. and Burkey L. A., 1956, J. Bacteriol., 71: 43-46.
3. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
4. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
5. Touraille C. and Bergere J. L., 1974, Biochimie, 56: 404-422.
6. Vissers M. M. M., Drichuis F., Giffel T., John P. D., Lankveld J. M. G., 2007, J. Dairy Sci., 90:928

Revision : 03/2020

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.