

Technical Data

Azotobacter Broth (Mannitol)

M1722

Azotobacter Broth (Mannitol) is recommended for cultivation of mannitol positive Azotobacter species from soil.

Composition**

Ingredients	Gms / Litre
Dipotassium phosphate	1.000
Magnesium sulphate	0.200
Sodium chloride	0.200
Ferrous sulphate	TRACE
Soil extract	5.000
Mannitol	20.000
Final pH (at 25°C)	8.3±0.2

^{**}Formula adjusted, standardized to suit performance parameters

Directions

Suspend 26.4 grams in 1000 ml distilled water. Heat if necessary to dissolve the medium completely. Dispense as desired. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Slight precipitate may occur after autoclaving, however it will not interfere with growth performance nor interfere with the interpretation of results.

Principle And Interpretation

Azotobacter is a free-living nitrogen-fixing bacterium (3), which is used as a biofertilizer in the cultivation of most crops. Azotobacter is found on neutral to alkaline soils, in aquatic environments, in the plant rhizosphere and phyllosphere. Azotobacter species are gram-negative aerobic soil-dwelling bacteria and are usually motile, oval, or spherical bacteria, form thick-walled cysts, and may produce large quantities of capsular slime. They are typically polymorphic, i.e. of different sizes and shapes. Their size of the cells ranges from 2-10 μm long and 1-2 μm wide. Plant needs nitrogen for its growth and Azotobacter fixes atmospheric nitrogen non-symbiotically. Therefore, all plants, trees, vegetables, get benefited. Beyond Azotobacter is used as a model it has biotechnological applications like use for alginate production and for nitrogen production in batch fermentations. This medium contains necessary nutrients for growth of Azotobacter species. For cultivation of mannitol positive Azotobacter species from soil, Azotobacter broth (Mannitol) can be used (1). It is used for cultivation of mannitol positive Azotobacter species from soil. It can also be useful for maintenance of Azotobacter species by adding extra 1% Mannitol to the medium containing agar i.e solid media as specified by the American Type Culture Collection (2).

Quality Control

Appearance

White to Cream homogeneous free flowing powder

Colour and Clarity of Prepared medium

Colourless to off-white coloured clear to slightly opalescent solution with slight precipitate forms intubes

Reaction

Reaction of 2.64% w/v aqueous solution at 25°C. pH: 8.3±0.2

pН

8.10-8.50

Cultural Response

Cultural characteristics observed after an incubation at 25-30°C for 24-48 hours.

Cultural Response

Organism	Inoculum (CFU)	Growth
Cultural Response		
Azotobacter beijerinckii	50-100	luxuriant
ATCC 12981		

HiMedia Laboratories Technical Data

Azotobacter nigricans ATCC 50-100 35009

luxuriant

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Reference

1.Pelczar M. Jr., 1957, Manual of Microbiological Methods.

2.ATCC Catalogue of Bacteria and Bacteriophages, 1992, 18th ed, American Type Culture Collection, Rockville, MD.

3.Subba Rao N. S., 1977, Soil Microorganisms and Plant Growth, Oxford and IBH Publishing Co., New Delhi.

Revision: 2 / 2015

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory,diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.