



Technical Data

Lysine Assay Medium

M1932

Intended Use:

Recommended for determining lysine concentration by microbiological assay method.

Composition**

Ingredients	Gms / Litre
Dextrose (Glucose)	50.000
Sodium acetate	40.000
Ammonium chloride	6.000
Potassium dihydrogen phosphate	1.200
Dipotassium hydrogen phosphate	1.200
Magnesium sulphate	0.400
Ferrous sulphate	0.020
Manganese sulphate	0.040
Sodium chloride	0.020
Adenine sulfate	0.020
Guanine hydrochloride	0.020
Uracil	0.020
Xanthine	0.020
Thiamine hydrochloride	0.001
Pyrodoxine hydrochloride	0.002
Pyridoxamine hydrochloride	0.600
Pyridoxal hydrochloride	0.600
Calcium panthothenate	0.001
Riboflavin	0.001
Nicotinic acid	0.002
p-Aminobenzoic acid	200.000mcg
Biotin	2.000mcg
Folic acid	20.000mcg
Glycine	0.200
DL-Alanine	0.400
Asparagine	0.800
L-Aspartic acid	0.200
L-Proline	0.200
DL-Serine	0.100
DL-Tryptophan	0.080
L-Glutamic acid	0.600
L-Histidine hydrochloride	0.124
DL-Phenylalanine	0.200
DL-Threonine	0.400
L-Tyrosine	0.200
DL-Valine	0.500
DL-Isoleucine	0.500
DL-Leucine	0.500
L-Arginine hydrochloride	0.484
L-Cystine	0.100
DL-Methionine	0.200
Final pH (at 25°C)	6.7±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 10.5 grams in 100 ml of purified / distilled water. Heat to boiling to dissolve the medium completely.

Mix well to distribute the slight precipitate evenly. Dispense in 5 ml amounts to each assay tube in increasing amounts of the

standard or the unknown and total volume 10 ml per tube is adjusted by addition of distilled / purified water. Sterilize by autoclaving at 15 lbs pressure (121°C) for 10 minutes.

Principle And Interpretation

Lysine Assay medium is formulated as described in Kavanagh (3). It contains all the essential growth factors for *Pediococcus acidilactici* ATCC 8042 except lysine. Exact concentration of lysine in the test material can be calculated by comparing results with standard curve of lysine.

Assay/Procedure: Stock cultures of the test organism *Pediococcus acidilactici* ATCC 8042, are prepared by stab inoculation into Lactobacilli Agar AOAC. Incubate the cultures at 35-37°C for 16-24 hours. After incubation, centrifuge aseptically and decant the supernatant. The pellet is washed 3-4 times with 10 ml of 0.85% NaCl solution. Then resuspend the cells in 10 ml of 0.85% NaCl solution. Inoculate each tube aseptically with 1 drop of the inoculum.

It is essential that a standard curve be set up for each assay since conditions of autoclaving, temperature of incubation, etc. which influence the standard curve readings, cannot be duplicated exactly from time to time.

Increasing amounts of the standard or the unknown and make up the volume to 10ml per tube containing 5 ml of the rehydrated medium. The growth response is measured turbidometrically.

Type of specimen

Isolated Microorganism

Specimen Collection and Handling

For isolated microorganism samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards. (3)

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

NA

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Off-white to yellow homogeneous free flowing powder

Colour and Clarity of prepared medium

Light amber coloured clear solution, which may contain a slight precipitate.

Reaction

Reaction of 10.5% w/v aqueous solution at 25°C. pH : 6.7±0.2

pH

6.50-6.90

Cultural Response

Microbiological Assay of Lysine was carried out using *Pediococcus acidilactici* ATCC 8042 after an incubation at 35-37°C for 16-20 hours .

Organism

Growth

Pediococcus acidilactici
ATCC 8042

Good growth is obtained. Gradual increase in growth with increasing conc.of standard L-Lysine 0, 30,60, 90, 120,150 mcg per assay tube was recorded as equivalent increase in absorbance at 660 nm.

Storage and Shelf Life

Store dehydrated and prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (1,2).

Reference

1. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
2. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
3. Kavanagh F., Analytical Microbiology Academic Press 1963, New York and London.

Revision : 01 / 2019

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.