

Technical Data

HiCrome[®] L. mono Differential Agar Base

Intended use

Recommended for the selective and differential isolation, enumeration and identification of *Listeria monocytogenes* and *Listeria* species based on PCPLC activity.

g / L 15.000 6.000 2.000 4.000 1.000 0.500 5.000 2.500 2.200 14.000

7.2±0.2

Composition**
Ingredients
Peptone
Tryptone
Yeast extract
Sodium pyruvate
Maltose
Magnesium glycerophosphate
Magnesium sulphate
Sodium chloride
Lithium chloride
Disodium hydrogen phosphate anhydrous
Chromogenic substrate
Agar

**Formula adjusted, standardized to suit performance parameters

Directions

Final pH (at 25°C)

Suspend 33.60 gram in 480 ml purified/distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Aseptically add sterile contents of 1 vial of Lecithin solution (FD332) and sterile rehydrated contents of PANC Selective Supplement (FD333). Mix well and pour into sterile Petri plates.

Principle And Interpretation

Listeria monocytogenes is a gram-positive foodborne human pathogen responsible for serious infections in pregnant women that may ultimately result in abortion, stillbirth, and birth of a child with neonatal listeriosis and meningitis or primary bacteremia in adults and juveniles (1). The pathogenicity of *Listeria ivanovii* for humans is uncertain. Since *L.monocytogenes* and *L.innocua* have similar biochemical properties, they cannot be differentiated on traditional media (PALCAM, Oxford). HiCrome[®] L.mono Differential Agar Base is based on for the selective and differential isolation of *Listeria* species on the basis of utilization of chromogenic substrate and lecithinase activity [Phosphotidylcholine phospholipase C (PCPLC)] (2). PI-PLC and PC-PLC, the major virulence factors, are only produced by pathogenic *L. monocytogenes* and *Listeria ivanovii* (3). Peptone, tryptone, yeast extract and sodium pyruvate provide nitrogenous and carbonaceous compounds, long chain amino acids, vitamins and essential growth nutrients . Maltose is the fermentable carbohydrate. Sodium chloride maintains osmotic equilibrium. Phosphate buffers the medium. Lithium chloride and added selective supplement (FD333) inhibit accompanying microflora and allow the growth of *Listeria* species. *Listeria* species hydrolyse the chromogenic substrate and produces green coloured colonies. Lecithin solution (FD332) helps in detecting PCPLC activity. Differentiation of *Listeria* species is based on phosphatidylcholine phospholipase C (PCPLC) activity. *L.monocytogenes* and *L.ivanovii* exhibits PCPLC activity which is seen as opaque halo around the colony.

Type of specimen

Food samples

Specimen Collection and Handling

For food samples, follow appropriate techniques for sample collection and processing as per guidelines (4). After use, contaminated materials must be sterilized by autoclaving before discarding.

M2009

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

- 1. Individual organisms differ in their growth requirement and may show variable growth patterns on the medium.
- 2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.
- 3. Further confirmation has to be carried using sensitivity testing.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.4% Agar gel

Colour and Clarity of prepared medium

Light amber coloured, opalescent gel forms in Petri plates

Reaction

Reaction of 6.72% w/v aqueous solution at 25°C. pH : 7.2±0.2

pН

7.00-7.40

Cultural Response

Cultural characteristics observed with added sterile PANC Selective Supplement (FD333) and Lecithin solution (FD332) after an incubation at 35 - 37°C for 24 - 48 hours.

Organism	Inoculum (CFU)	Growth	Recovery	Colour of Colony	PCPLC activity
Escherichia coli ATCC 25922 (00013*)	>=10 ⁴	inhibited	0%		
Escherichia coli ATCC 8739 (00012*)	>=10 ⁴	inhibited	0%		
Enterococcus faecalisATCC 29212 (00087*)	£ >=10 ⁴	inhibited	0%		
Listeria innocua ATCC 33090 (00017*)	50-100	luxuriant	>=50%	greenish-blue	negative
Listeria ivanovii subsp. ivanovii serovar 5 ATCC 19119 (00018*)	50-100	luxuriant	>=50%	greenish-blue	positive, opaque halo around the colony exhibiting phophatidylcholine phospholipase acivity
Listeria monocytogenes serovar l ATCC 19112	50-100	luxuriant	>=50%	greenish-blue	positive, opaque halo around the colony exhibiting phosphatidylcholine specific phospholipase activity
Listeria monocytogenes serovar 4b ATCC 13932 (00021*)	50-100	luxuriant	>=50%	greenish-blue	positive, opaque halo around the colony exhibiting phosphatidylcholine specific phospholipase activity
Listeria monocytogenes serovar l/2a ATCC 35152 (00109*)	50-100	luxuriant	>=50%	greenish-blue around	positive, opaque halo the colony exhibiting phosphatidylcholine specific phospholipase activity

Key : (*) - Corresponding WDCM Numbers

Storage and Shelf Life

Store between 15-25°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (5,6).

Reference

- 1. Painter J, Slutsker L. 2007. Listeriosis in humans, p 85–109. In Ryser ET, Marth EH (ed), *Listeria*, listeriosis, and food safety. Marcel Dekker, New York, NY.
- 2. Sang-Hyun Park, Pahn-Shick Chang, Sangryeol Ryu and Dong-Hyun Kang. Development of a Novel Selective and Differential Medium for the Isolation of *Listeria monocytogenes*. Applied and Environmental Microbiology 2014.
- Mengaud J, Braun-Breton C, Cossart P 1991. Identification of phosphatidylinositol-specific phospholipase C activity in Listeria monocytogenes: a novel type of virulence factor. Mol. Microbiol. 5:367–372. doi:10.1111/j.1365-2958.1991
- 4. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
- 5. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
- Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Revision: 05/2024

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.