

Technical Data

Sugar Free Agar

M307

Sugar Free Agar is recommended for the examination of butter in accordance with International Dairy Federation.

Composition**	
Ingredients	Gms / Litre
Casein enzymic hydrolysate	7.500
Pancreatic digest of gelatin	7.500
Sodium chloride	5.000
Agar	14.000
Final pH (at 25°C)	7.6±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 34 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and pour into sterile Petri plates.

Principle And Interpretation

Butter may be contaminated in the manufacturing process, from contaminated water sources, air or unsanitary equipment. Yeast, moulds and coliforms do not survive pasteurization and are rarely found in butter made under sanitary conditions. Estimates of the number of coliforms, yeasts and moulds or Enterococci in samples of butter taken at various stages of processing are useful in tracing the source of contamination (1). Sugar Free Agar is prepared as recommended by the International Dairy Federation (2) and is based on the formulation developed by Ritter and Eschmann (3) for the detection and enumeration of contaminating bacteria in butter and other processed dairy products.

The medium contains pancreatic digest of gelatin and casein enzymic hydrolysate which provide essential nutrients to the microorganisms. The medium is free of any fermentable carbohydrate and has relatively little nutrient value. Therefore microorganisms that do not form the normal flora of butter and other processed milk products can be grown selectively on this medium. Psychrotrophic and mesophilic gram-negative bacteria also grow well on this medium. Some Lancefield group D streptococci are not able to grow on this medium and hence do not interfere in the total count of the gram-negative bacteria (3). This medium aids differentiating non-lactic acid contaminants from lactic acid bacteria.

Quality Control

Appearance Cream to yellow homogeneous free flowing powder Gelling Firm, comparable with 1.4% Agar gel Colour and Clarity of prepared medium Medium amber coloured clear to slightly opalescent gel forms in Petri plates Reaction Reaction of 3.4% w/v aqueous solution at 25°C. pH : 7.6±0.2 pH 7.40-7.80 **Cultural Response** M307: Cultural characteristics observed after an incubation at 35-37°C for 18-48 hours. Organism Inoculum Growth Recovery (CFU)

Cultural Response

Please refer disclaimer Overleaf.

Escherichia coli ATCC 25922	50-100	good-luxuriant >=70%
Bacillus cereus ATCC 10876 Candida albicans ATCC	50-100 50-100	good-luxuriant >=70% good-luxuriant >=70%
10231 Enterococcus faecalis ATCC	50-100	good-luxuriant >=70%
29212 Leuconostoc mesenteroides	50-100	good-luxuriant >=50%
ATCC 8293		-

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Reference

1. International Dairy Federation, 1964, International Standards FIL-IDF 30

2. Ritter P. and Eschmann K. H., 1966, Alimenta, 5 (2): 433. Thomas S. B., 1969, J. Appl. Bacteriol., 32: 269

Revision : 02 / 2015

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516, Swastik Disha Business Park, Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com