



Algae Culture Broth

M342

Intended Use:

Recommended for cultivation of algae from soil, water and sewage. Also for preparing the inoculum for the bioassay of algicidal chemicals.

Composition**

Ingredients	Gms / Litre
Sodium nitrate	1.000
Dipotassium hydrogen phosphate	0.250
Magnesium sulphate	0.513
Ammonium chloride	0.050
Calcium chloride	0.058
Ferric chloride	0.003
Final pH (at 25°C)	7.0±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 1.87 grams in 1000 ml purified / distilled water. Heat if necessary to dissolve the medium completely. Dispense in tubes or flasks as desired. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

Principle And Interpretation

Algae (singular alga) encompass several groups of relatively simple living aquatic organisms that capture light energy through photosynthesis, using it to convert inorganic substances into organic matter. Algae range from single-cell organisms to multicellular organisms, some with fairly complex differentiated form and (if marine) called seaweeds. Algae are usually found in damp places or water bodies and thus are common in terrestrial as well as aquatic environments. Various algae play significant roles in aquatic ecology. Algae are used by humans in a number of ways. Because many species are aquatic and microscopic, they are cultured in clear tanks or ponds and either harvested or used to treat effluents pumped through ponds (2,5). Algae Culture Broth is recommended for the isolation and cultivation of algae from soil, water and sewage. Algae Culture Broth is used to prepare the inoculum for the bioassay of algicidal chemicals. Algae Culture Broth is similar in composition to Algae Culture Agar, except the agar.

The medium provides all necessary nutrients for good growth of Algae but does not provide for other than minimal growth of bacteria and fungi.

Type of specimen

Water samples: Wastewater, Soil sample

Specimen Collection and Handling

For water samples follow appropriate techniques for handling specimens as per established guidelines (3,4). For soil samples, follow appropriate techniques for sample collection and processing as per guidelines (6). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

1.This medium is general purpose medium and may not support the growth of fastidious organisms.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

White to light yellow homogeneous free flowing powder

Colour and Clarity of prepared medium

White coloured clear to slightly opalescent solution in tubes.

Reaction

Reaction of 0.187% w/v aqueous solution at 25°C. pH : 7.0±0.2

pH

6.80-7.20

Cultural Response

Cultural characteristics observed under suitable light source after an incubation at 20-25°C within 1 week.

Organism

Growth

Chlorella pyrenoidosa good-luxuriant
ATCC 50476

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 15-25°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (3,4).

Reference

1. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C.
2. Guiry M. D. and Blunden G., (Ed.), 1991, Seaweed Resources in Europe: Uses and Potential. John Wiley and Sons Ltd.
3. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
4. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
5. Lembi C. A. and Waaland J. R., (Ed.), Algae and Human Affairs, 1988, Cambridge University Press.
6. Subba Rao N. S., 1977, Soil Microorganisms and Plant Growth, Oxford and IBH Publishing Co., New Delhi.

Revision :02 / 2020

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.