



Buffered Yeast Agar

M585

Intended Use:

Recommended for cultivation of yeasts and moulds and for the controlling of bottle washing operations in the soft drinks and related industries.

Composition**

Ingredients	Gms / Litre
Yeast extract	5.000
Dextrose (Glucose)	20.000
Ammonium sulphate	0.720
Ammonium dihydrogen phosphate	0.260
Agar	15.000
Final pH (at 25°C)	5.5±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 41 grams in 1000 ml purified / distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 115°C for 20 minutes. Mix well and pour into sterile Petri plates.

Principle And Interpretation

Yeasts grow well on a minimal medium containing only dextrose and salts. The addition of yeast extract allows faster growth so that during exponential or log phase growth, the cells divide every 90 minutes (1). Buffered Yeast Agar is prepared as per the modification of the yeast-salt medium described by Davis (4).

The medium contains yeast extract, which supplies B-complex vitamins to stimulate growth. Dextrose is the carbohydrate source. The reaction of this medium can be adjusted to required pH values by the addition of citric or lactic acid to the medium after sterilization. The following table shows the amount of the acids required to be added to 100 ml of Buffered Yeast Agar cooled to 50°C.

Volume of acid to be added to 100 ml of medium to achieve the desired pH

pН	1% w/v solution of	1% w/v solution of	
	Citric acid monohydrate (ml)	Lactic acid (ml)	
4.75	1.26	0.125	
4.5	2.24	0.2	
4.25	3.92	0.3	
4.0	6.16	0.45	
3.75	9.52	0.7	
3.5	14.56	1.17	

Type of specimen

Bottles used in soft drink industries.

Specimen Collection and Handling:

Bunker (2,3) described a practical method for assessing the efficiency of the bottle cleaning operations. In this method, the bottle under test is converted into a roll-tube culture by coating it internally with the medium. When the agar sets, the bottle is incubated and the colonies are counted and examined. This method gives better results than rinsing the bottle and subsequently plating the rinsings. When used for this purpose, the agar concentration in Buffered Yeast Agar should be increased by 1%w/v (before sterilization).

Warning and Precautions :

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per

established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations :

1. This medium is general purpose medium and may not support the growth of fastidious organisms.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Light amber coloured, clear to slightly opalescent gel forms in Petri plates.

Reaction

Reaction of 4.1% w/v aqueous solution at 25°C. pH : 5.5±0.2

pН

5.30-5.70

Cultural Response

Cultural characteristics observed after an incubation at 25-30°C for 48-72 hours.

Organism	Inoculum (CFU)	Growth	Recovery
Candida albicans ATCC 10231 (00054*)	50-100	good-luxuriant	>=70%
Saccharomyces cerevisiae ATCC 9763 (00058*)	50-100	good-luxuriant	>=70%
#Aspergillus brasiliensis ATCC 16404 (00053*)	50-100	good-luxuriant	

Key:- *Corresponding WDCM numbers. #Formerly known as *Aspergillus niger*

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 20-30°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (5,6).

Reference

1. Ausubel, Brent, Kingston, Moore, Seidman, Smith and Struhl, 1994, Current Protocols in Molecular Biology, Current Protocols, Brooklyn, N.Y.

- 2. Bunker H. J., 1952, Lab. Prac., 18:354.
- 3. Bunker H. J., 1956, Wallerstein Lab. Communications, 19(65): 143.
- 4. Davis J. G., 1931, J. Dairy Res., 3:133.
- 5. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
- Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg.office : 23, Vadhani Ind.Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6116 9797 Corporate office : A-516,Swastik Disha Business Park,Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com