

Technical Data

Lactobacillus MRS Agar(MRS Agar)

Recommended for the isolation and enumeration of lactic acid bacteria from meat and meat products. The composition and performance criteria of this medium are as per the specifications laid down in ISO 13721:1995.

Composition**	
Ingredients	Gms / Litre
Meat extract B [#]	8.000
Peptone	10.000
Yeast extract	5.000
Ammonium citrate	2.000
Sodium acetate	5.000
Magnesium sulphate, heptahydrate	0.200
Manganese sulphate, tetrahydrate	0.050
Dipotassium phosphate	2.000
Glucose, anhydrous	20.000
Polysorbate 80	1.000
Agar	12.000
Final pH (at 25°C)	5.7±0.2
**Formula adjusted, standardized to suit performance parameters	

- Equivalent to Beef extract

Directions

Suspend 65.13 grams (the equivalent weight of dehydrated medium per litre) in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle And Interpretation

Lactobacilli MRS medium is based on the formulation of deMan, Rogosa and Sharpe (1) with slight modification. It supports luxuriant growth of all Lactobacilli from oral cavity (1), dairy products (2), foods (3), faeces (4) and other sources (5). Lactobacillus MRS Agar is recommended by ISO Committee (6).

Peptone and meat extract B supplies nitrogenous and carbonaceous compounds, long chain amino acids and other essential growth nutrients. Yeast extract provides vitamin B complex and glucose is the fermentable carbohydrate and energy source. Polysorbate 80 supplies fatty acids required for the metabolism of Lactobacilli. Sodium acetate and ammonium citrate inhibit Streptococci, moulds and many other microorganisms. Phosphates provide good buffering action in the media.

Lactobacilli are microaerophillic and generally require layer plates for aerobic cultivation on solid media. When the medium is set, another layer of un-inoculated MRS Agar is poured over the surface to produce a layer plate (5). Lactobacilli isolated on MRS Agar should be further confirmed biochemically.

Type of specimen

Food and dairy samples ; Water samples

Specimen Collection and Handling:

For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (6,7). For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (1). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

Further bichemical tests are required for confirmation.

M641I

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to light yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.2% Agar gel.

Colour and Clarity of prepared medium

Medium to dark amber coloured slightly opalescent gel forms in Petri plates

Reaction

Reaction of 6.51% w/v aqueous solution at 25°C. pH : 5.7±0.2

pН

```
5.50-5.90
```

Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours (longer if neccesary).(with 5% CO2) **Cultural Response**

Organism	Inoculum (CFU)	Growth	Recovery
Cultural Response Lactobacillus acidophilus ATCC 4356	50-100	luxuriant	>=50%
Lactobacillus casei ATCC 9595	50-100	luxuriant	>=50%
Lactobacillus fermentum ATCC 9338	50-100	luxuriant	>=50%
Lactobacillus plantarum ATCC 8014	50-100	luxuriant	>=50%
Lactococcus lactis subsp. lactis ATCC 19435	50-100	luxuriant	>=50%
Lactococcus sakei ATCC 15521	50-100	luxuriant	>=50%
Pediococcus damnosus ATCC 29358	50-100	luxuriant	>=50%
Pediococcus pentosaceus ATCC 33316	50-100	luxuriant	<=10%
<i>Bifidobacterium bifidum*</i> ATCC 11863	50-100	luxuriant	<=10%
Escherichia coli ATCC 25922	50-100	none-poor	<=10%

Key: * -Growth under anaerobic conditions for 72 hours

Storage and Shelf Life

Store below 8°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Reference

1.deMan J., Rogosa M. and Sharpe M., 1960, J. Appl. Bacteriol., 23:130.

2. Marshall R.T. (Ed.), 1992, Standard Methods for the Examination of Dairy Products, 16th ed., APHA, Washington, D.C.

3.Downes F. P. and Ito K., (Eds.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 4th Ed., APHA, Washington, D.C.

4.Sabine and Vaselekos, 1965, Nature, 206:960.

5.MacFaddin J., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol.1, Williams and Wilkins, Baltimore.

6. International Organization for Standardization (ISO), 1995, Draft ISO/DIS, 13720.

Revision : 05/2023

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Corporate Office : Plot No.C-40, Road No.21Y, MIDC, Wagle Industrial Area, Thane (W) - 400604, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com