
Tween® 20, For Molecular Biology

MB067

Product Identifier

CAS No.	:	9005-64-5
EC No.	:	500-018-3
Synonym	:	Polyoxyethylenesorbitan monolaurate; Polysorbate-20
HS Code	:	3402 42 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Lemon to amber liquid
Solubility	:	33.3 mg soluble in 1 mL of water
FTIR	:	Matches with the standard pattern
DNases	:	None detected
RNases	:	None detected
Hydroxyl value	:	96 - 108
Saponification value	:	40 - 50
Acid value	:	<= 2.0
Residue on ignition	:	<= 0.25%
Water (K.F.)	:	<= 3.0%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	TR7400000
WGK	:	1