
DNA sodium salt, From calf thymus, highly polymerised, For Molecular Biology

MB103

Product Identifier

| | | |
|------------|---|-----------------------------------|
| CAS No. | : | 73049-39-5 |
| Synonym | : | Deoxyribonucleic acid sodium salt |
| HS Code | : | 2934 99 90 |
| Storage | : | 2 - 8°C |
| Shelf life | : | 3 years |

Technical Specification

| | | |
|------------|---|--|
| Appearance | : | White to off-white to brown powder or threads or solid |
| Solubility | : | 10 mg soluble in 10 mL of water |
| DNases | : | None detected |
| RNases | : | None detected |

Safety Information

| | | |
|---------------|---|---------------------|
| UN No. | : | Not dangerous goods |
| Class | : | - |
| Packing Group | : | - |
| RTECS | : | Not available |
| WGK | : | 1 |