

Sarcosyl, For Molecular Biology

MB108


Product Identifier

CAS No.	:	137-16-6
EC No.	:	205-281-5
Molecular Formula	:	C ₁₅ H ₂₈ NO ₃ Na
Molecular Weight	:	293.38
HS Code	:	2924 19 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to beige powder
Solubility	:	33.3 mg soluble in 1 mL of water
DNases	:	None detected
RNases	:	None detected
FTIR	:	Matches with the standard pattern
Water (K.F.)	:	<= 7.00%
Assay (HPLC/NT, anhydrous basis)	:	94.00 - 106.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H315- H318- H330
Precautionary Statement(s)	:	P260- P280- P304 + P340- P305+P351+P338- P310- P403 + P233
UN No.	:	2811
Class	:	6.1
Packing Group	:	II
RTECS	:	Not available
WGK	:	1