
DEAE Cellulose - 52, For Molecular Biology

MB110

Product Identifier

CAS No.	:	9013-34-7
Synonym	:	(Diethylaminoethyl-cellulose) fine mesh For column chromatography
HS Code	:	3914 00 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to off-white powder or solid
Solubility	:	Insoluble in water
DNases	:	None detected
RNases	:	None detected
Water (K.F.)	:	min. 65.00%
Capacity	:	0.80 - 1.10 meq/g, dry weight

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1